



ANALYTICAL INSTRUMENT and SMALL EQUIPMENT MANAGEMENT AND CONTROL

Each piece of equipment or instrument is barcoded and/or RFID tagged. Location (lab room number) can be associated with barcode and/or RFID tag.



A RFID reader is placed in every room of the laboratory



Agencies requiring to have an inventory of equipment by location, can accomplish this with a push of a button

Agencies can find the whereabouts of small pieces of equipment that are due for service or calibration (such as scales and pipettes).

By having a RFID reader at the main entrance/exit of lab, laptops entering or leaving the laboratory is recorded.